


Integrated Clinical Environment (ICE) Overview:

The ICE system contains simulated devices that monitor patient data, a simulated infusion pump that monitors a drug to a patient, and a supervisor app that displays patient vitals, alarms, and is used to control the valid ranges for the system. For example, the supervisor app can decide what a valid pulse rate is for a patient. The supervisor app and the device apps can both show alarms in certain circumstances.



What's in the system?

Simulated devices:

- Simulated devices send their device ID, an image representing the device, and numeric data such as patient vital signs. Simulated devices receive alarm limits that tell them the ranges of vital sign values that should produce an alarm (in the case where the device can display an alarm).

Devices send and receive much more than this, but for simplicity we will describe only some of the important values.

Infusion pump:

- The infusion pump produces and consumes a lot of the same data as the other devices, but it also needs to be able to receive a command telling it to stop infusing immediately.


Supervisor:

- The supervisor receives the device IDs, images, and vital signs from the devices. It receives the status of the infusion pump.
- The supervisor sends and receives alarm limits, which are used to define what the valid ranges are for each vital sign. Both the supervisor and some devices use these values to decide whether to display an alarm. Or not. This is what we will compromise.

Exercise 1:


Overview:

The ICE demo allows you to create simulated devices, and to view them and control their behavior. In this exercise, you will run the ICE demo with no security enabled at all. The goal of this exercise is to understand the data in the ICE system, to visualize it over the network, and to understand what a “sniffer” app (legitimate or otherwise) can see.



Exercise:

1. Open Wireshark, and start capturing traffic. Double-click on the “Ethernet” button and you will start to see packets. (If you don’t see any traffic right away, you will see some when you start the ICE application.)




Capturing from Ethernet (not tcp port 3389)						
File Edit View Go Capture Analyze Statistics Telephony Help						
Apply a display filter ... <Ctrl-/>						
No.	Time	Source	Destination	Protocol	Length	Info
1	0.000000	Vmware_b4:60:d7	Vmware_30:34:a0	ARP	60	Who has 10.160.196.144? Tell 10.160.0.1
2	0.000021	Vmware_30:34:a0	Vmware_b4:60:d7	ARP	42	10.160.196.144 is at 00:50:56:30:34:a0
3	8.686084	fe80::ffff:ffff:ffff:ff02::2	ff02::2	ICMPv6	103	Router Solicitation
4	8.799672	fe80::8000:f227:a10...	fe80::ffff:ffff:ffff:ffff	ICMPv6	151	Router Advertisement
5	16.047844	10.160.196.144	10.160.255.255	NBNS	92	Name query NB WPAD<00>
6	16.048077	10.160.196.144	224.0.0.252	LLMNR	64	Standard query 0x0225 A wpad
7	16.048232	10.160.196.144	224.0.0.252	LLMNR	64	Standard query 0x3cec AAAA wpad

▶ Frame 1: 60 bytes on wire (480 bits), 60 bytes captured (480 bits) on interface 0
▶ Ethernet II, Src: Vmware_b4:60:d7 (00:50:56:b4:60:d7), Dst: Vmware_30:34:a0 (00:50:56:30:34:a0)
▶ Address Resolution Protocol (request)


```
0000 00 50 56 30 34 a0 00 50 56 b4 60 d7 08 06 00 01 .PV04..P V.`.....
0010 08 00 06 04 00 01 00 50 56 b4 60 d7 0a a0 00 01 .....P V.`.....
0020 00 00 00 00 00 00 0a a0 c4 90 00 00 00 00 00 .....
0030 00 00 00 00 00 00 00 00 00 00 00 00 ..... .....
```

You may want to filter for RTPS traffic if you would like to see only the data specific to our application:




2. Go to the OpenICE-0.6.3\bin. Run the OpenICE.bat file.

(**Note:** for the purposes of these exercises, you must always run the open ICE project from within the <exercise>\OpenICE-0.6.3\bin directory, either by double-clicking it from the windows GUI or changing into that directory at the command-prompt.)





3. Keep the application in domain 0 (the default)

4. Click on "Start ICE_Supervisor"



5. Once the full application comes up, click on "Create an ICE Device Adapter" in the lower right hand corner and create a simulated pulse oximeter. (This has easy data to view).



6. Run the c:\Users\student\ndds.5.1.1\scripts\nddsspy tool on another machine,

with the option -printSample. This tool is a debugging / simple visualization tool. Without encryption, you can visualize all the data in the system. By running rtiddsspy on a different machine, you will be able to see all the traffic in Wireshark.

7. View the RTPS traffic using Wireshark. This traffic includes discovery data that describes the data being written and read by each application.

What you should see:

Packets that are marked in blue in Wireshark are discovery data, and they show the metadata about what is being done in the system.

14	5.369795	10.160.196.1	239.255.0.1	RTPS	702 INFO_TS, DATA(p)
16	5.395194	10.160.196.1	239.255.0.1	RTPS	1418 INFO_TS, DATA(r)
18	5.397431	10.160.196.1	239.255.0.1	RTPS	1342 INFO_TS, DATA(r)
19	5.458231	10.160.196.1	239.255.0.1	RTPS	1442 INFO_TS, DATA(r)

Packets that are marked in red are the runtime data of the system.

464	44.269828	10.160.196.1	10.160.86.31	RTPS	158 INFO_TS, DATA
465	44.270214	10.160.196.1	10.160.86.31	RTPS	174 INFO_TS, DATA
466	44.295673	10.160.196.1	10.160.86.31	RTPS	110 INFO_DST, HEARTBEAT
467	44.295910	10.160.196.1	239.255.0.1	RTPS	146 INFO_TS, DATA(m)

If you have created a pulse oximeter, some of its data contains the string MDC_PULS_OXIM_PULS_RATE (describing the pulse rate data), and the string MDC_DIM_DIMLESS (saying this data does not have dimensions). You can search for one of these values in Wireshark to see the data as it is being sent.

The user data transferred in a ISSUE submessage (rtsp.issueData), 96 bytes